

B¹
adult cells, as donor nuclei in the nuclear transfer procedure. The dogma has been that embryonic cells are more easily reprogrammed than adult cells. This dates back to earlier NT studies in the frog (review by DiBerardino, *Differentiation*, 17:17-30 (1980)). Wolfe et al. also did not report that the inner cell mass cells from NT units could be used to form an ES cell-like colony that could be propagated. - -

Please replace the paragraph beginning at line 22 of page 19 with the following amended paragraph:

B²
- - The activated NT units then may be cultured in a suitable *in vitro* culture medium until the generation of embryonic or stem-like cells and cell colonies. Culture media suitable for culturing and maturation of embryos are well known in the art. Examples of known media, which may be used for bovine embryo culture and maintenance, include Ham's F-10 + 10% fetal calf serum (FCS), Tissue Culture Medium-199 (TCM-199) + 10% fetal calf serum, Tyrodes-Albumin-Lactate-Pyruvate (TALP), Dulbecco's Phosphate Buffered Saline (PBS), Eagle's and Whitten's media. One of the most common media used for the collection and maturation of oocytes is TCM-199, and 1 to 20% serum supplement including fetal calf serum, newborn serum, estrual cow serum, lamb serum or steer serum. A preferred maintenance medium includes TCM-199 with Earl salts, 10% fetal calf serum, 0.2 mM Na pyruvate and 50 µg/ml gentamicin sulphate. Any of the above may also involve co-culture with a variety of cell types such as granulosa cells, oviduct cells, BRL cells, uterine cells and STO cells. - -
